Primospine

Ortesis sacrolumbar y toraco-lumbosacro rígida
 جهاز التقويم الصلب للفقرات القطنية العجزية والصدروالفقرات القطنية العجزية
 Pevná lumbosakrální a torakolumbosakrální ortéza
 Starre sakrolumbale und thorakolumbosakrale Orthese
 Rigide lumbosakral og thorako-lumbosakral ortoser
 Άκαμπτος κηδεμόνας οσφύος και θωρακοοσφυϊκός
 Rigid thoracolumbosacral and sacrolumbar orthosis
 Jäykkä ristiselkäortoosi ja jäykkä rintaranka-ristiselkäortoosi
 Orthèse sacro-lombaire et thoraco-lombo-sacré rigide
 Kruta lumbosakralna i torakolumbosakralna ortoza
 Merev mellkas-ágyék-keresztcsont ortézis
 Ortesi rigida lombosacrale e toracico-lombo-sacrale
 Stijve sacro-lumbale en thorax-lumbosacro orthese
 Stiv lumbosakral og thorax-lumbosakral ortose
 Szttywne ortczy lędźwiowo-krzyżowa i piersiowo-lędźwiowo-krzyżowa
 Ortótese sacrolombar e toracossacrolombar rígida
 Orteză sacro-lombară și toraco-lombo-sacrală rigidă
 Жесткий пояснично-крестцовый и грудно-пояснично-крестцовый ортез
 Pevná lumbosakrálna a torako-lumbosakrálna ortéza
 Rigid ländryggs- och thoraco-lumbosakralortos
 Trda lumbosakralna in torako-lumbosakralna ortoza

PRS660 · PRS661
 PRS666

Instrucciones de adaptación · تعليمات المستخدم · Pokyny K Nastavení · Anleitung ·
 Tilpasning instruksjon · Οδηγίες εφαρμογής · Fitting instructions · Sovitusohjeet ·
 Instructions d'adaptation · Upute za postavljanje · Rögzítési utasítások · Istruzioni
 per l'adattamento · Montage-instructies · Tilpassingsinstruksjon · Instrukcja obsługi
 · Instruções de adaptação · Instructiuni de reglare · Инструкции по адаптации ·
 Návod na použitie · Utprovningsinstruktion · Navodila za prilagoditev

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Description and intended use:

Rigid sacrolumbar orthosis (PRS660-PRS661): The rigid sections are made of high-density and low-density thermoplastics, with highly breathable and comfortable materials that provide the patient with a high level of sacrolumbar stabilisation.

Rigid thoracolumbosacral orthosis (PRS666): The rigid sections are made of high-density and low-density thermoplastics, with highly breathable and comfortable materials that provide the patient with a high level of thoracolumbar stabilisation.

Indications:

Rigid sacrolumbar orthosis: spinal disc herniation, osteoporosis, spondyloarthritis, spondylolysis (with or without stability), trauma, compression fractures, degenerative disc disease, postoperative fusion and chronic muscle disease.

Rigid thoracolumbosacral orthosis: spinal disc herniation, osteoporosis, spondyloarthritis, spondylolysis (with or without stability), trauma, compression fractures, degenerative disc disease, postoperative fusion, chronic muscle disease, and prevention of canal stenosis in bedridden patients.

Contraindications: Not indicated for use on open wounds with swelling, reddening, or heat buildup.

Warnings:

- Read these instructions carefully before using the orthosis. If you have any questions, consult your doctor or the shop that sold you the product.
- For the orthosis to be effective, patients must use the appropriate size. If the patient is close to the maximum range of a size and has a protruding abdomen (but not a pendulous abdomen), choose the next size up.
- Always follow these general instructions of use and any particular guidelines indicated by your general practitioner.
- Only a doctor is qualified to prescribe and determine the duration of treatment and carry out subsequent monitoring.
- If you notice any side effects, skin conditions or sensitivity, consult your doctor immediately. Discontinue use until assessed by a specialist.
- Your doctor might consider the use of topical creams in conjunction with the orthosis.
- In areas where pressure is applied, the skin should not be damaged or hypersensitive.
- In the event of mild irritation caused by sweat, we recommend wearing a cotton t-shirt to act as a buffer between the skin and the fabric.
- Although the orthosis is not for single use, it should only be used by one patient. The product should not be reused by other patients.
- The first fitting and adjustment must be carried out by a qualified professional (doctor, orthopaedic technician, etc.).
- For the orthosis to function correctly, its components must be regularly checked to ensure they are in good condition. The medical staff responsible for monitoring the treatment may recommend the product to the patient or suggest that it be replaced if it is damaged or worn.
- To ensure they function correctly, clean the hook-and-loop fastenings regularly, removing any material that may have stuck to them.
- Some of the materials that make up the product are not flame retardant. Therefore, its use is not recommended in environments with excessive heat, fire or radiation.
- Clean the product regularly and maintain proper hygiene.
- Ask your doctor if you can remove the orthosis during bathing or showering. If you cannot, make sure you cover it or keep it away from the water.
- Do not let children play with this device.

Instructions for fitting (PRS660-PRS661):

1. Position the rear section of the orthosis on the patient's back at the height of the sacrolumbar area, making sure it is centred as shown in fig. 1. Leave the side reinforcement straps loose, so that they do not impede the adjustment of the orthosis. Then centre the front section of the orthosis over the patient's abdomen.
2. If necessary, and before closing the orthosis, shape the stays to fit the anatomical curves of the patient.
3. Close the orthosis by fastening the microhook side sections to the front area, so that the orthosis is completely adapted to the patient's body shape as shown in fig. 2.
4. Lastly, adjust the side reinforcement straps in order to achieve the desired level of compression, then fasten the straps to the front of the orthosis. Adjust the lower straps to the required tension

and then the upper ones, until the orthosis is perfectly adapted to the patient as shown in fig. 3.

Instructions for fitting (PRS666):

1. Position the rear section of the orthosis on the patient's back, followed by the retropulsion straps. Make sure that the rear section is centred. To position the straps correctly, pass them beneath the patient's underarms, then up over the shoulders to the back again. Then insert the straps into the fastener located in the upper-rear section of the back part of the orthosis, as shown in fig. 4.
2. Leaving all of the straps loose, so that they do not impede the adjustment of the orthosis, centre the front section of the orthosis over the patient's abdomen.
3. If necessary, and before closing the orthosis, shape the stays to fit the anatomical curves of the patient.
4. Close the orthosis by fastening the microhook side sections to the front area, so that the orthosis is completely adapted to the patient's body shape as shown in fig. 5.
5. Adjust the side reinforcement straps in order to achieve the desired level of compression, then attach the straps to the front of the orthosis. Adjust the lower straps to the required tension and then the upper ones, until the orthosis is perfectly adapted to the patient as shown in fig. 6.
6. Lastly, adjust the retropulsion straps, making sure that they cross each other at the back to form an X shape, as shown in fig. 4. Then slide them towards the sides, in the direction of the anterior abdominal wall. Tension the straps with both hands and attach them to the front of the orthosis using the microhooks, as shown in fig. 7, until the desired level of retropulsion is achieved.

Note: The straps are fitted with crocodile-type microhook fasteners at the ends, so that - where necessary - the length of the straps can be adapted to the patient's height and the positioning of the orthosis during the first fitting and adjustment. If you detach the fasteners, the end of the strap will be left free and the excess length can be trimmed. Then reattach the fasteners onto both straps.

Adapting the thermoplastic inserts:

The orthosis is supplied with two thermoplastic inserts: a high-density lumbar insert and a low-density abdominal insert. The abdominal insert does not require adaptation. However, if you need to adapt the lumbar insert, follow these instructions:

1. Remove the thermoplastic insert from its fabric cover.
2. Adapt the insert to the patient's anatomy using a heat gun at 180°C, maintaining a distance of 10-15 cm.
3. Once the insert has been adapted, place it back inside the fabric cover.

Instructions for adapting the orthosis to patients with a pendulous abdomen:

For patients with a pendulous abdomen, the abdominal wall must be lifted after the orthosis has been fitted. To do this, the lower strap must be attached as low as possible above the pubis, and fastened first. The upper strap can then be tensioned and attached. If the model only has a single strap, this must be attached as low as possible above the pubis.

Product maintenance - washing instructions:

Thermoplastic inserts: Remove the plastic inserts from the fabric covers. They can be cleaned with a damp cloth.

Orthosis: *Before washing, remove the thermoplastic valves and make sure all of the microhook fasteners are closed. Wash with a neutral detergent at a maximum of 40°C. Do not use bleach. Do not tumble dry. Do not iron. Do not dry clean. Do not expose to direct sources of heat such as stoves, heaters, radiators, direct sunlight, etc.

Note: Any serious incident involving the product must be reported to the manufacturer and to the competent authority of the Member State where the user and/or patient is located.

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PRIM, S.A.